



LabCorp Employer Services COVID-19 Pixel by LabCorp™ Testing Experience

has partnered with LabCorp Employer Services (LES) to provide COVID-19 testing to support a healthy and safe work environment for all employees. The Pixel by LabCorp™ at-home COVID-19 test detects the presence of the SARS CoV-2 virus via a self-collected nasal swab. In an effort to support you with at-home testing, we've outlined what you can expect when ordering and completing this process.



To order your test kit, use the LES secure online HIPAA compliant portal, WellConnect Plus:



Kit shipped overnight

You will be directed to the Pixel by LabCorp website where you can provide shipping information in order to receive your kit. You will be able to track your kit within the Pixel site. If you order your test kit before 3 PM EST, it will be sent overnight via FedEx and arrive the next day. If you order after 3 PM EST, it will be shipped the following morning.



The kit includes all the supplies you will need, including detailed instructions to complete your self-collection.



Following the instructions included in the kit, you will take the provided cotton swab out of the package, taking care to not touch the cotton tip. You will then need to rotate the swab and insert less than one inch into each nostril, rotating at least three times per nostril.



Once you have completed the self-collection, you will place your sample in the prepaid, pre-addressed FedEx return package and deliver it to a FedEx drop box, per the guidelines in the test kit instructions. You will be able to track your sample within the Pixel by LabCorp site.



Your results will be available in two to four days on the secure online, HIPAA compliant portal where you first requested your kit. Should your results indicate positive for COVID-19, you will be contacted by a health care professional to help you with next steps.



www.LabCorp.com